## PHOSPHORUS-SL ASSAY

**CATALOGUE NUMBER:** 117-30 **SIZE:** R1: 3 x 100 mL, R2: 1 x 75 mL

### INTENDED USE


For the IN VITRO quantitative measurement of inorganic phosphorus in serum.

### TEST SUMMARY

Increased serum inorganic phosphorus concentrations (hyperphosphatemia) are usually the result of Vitamin D overdose, hypoparathyroidism, or renal failure. Decreased serum concentrations (hypophosphotemia) usually result from rickets, hyperparathyroidism, or Fanconi Syndrome (a defect in absorption of phosphorus and other metabolites from the glomerular filtrate).<sup>(1)</sup>

The majority of methods for determining inorganic phosphorus are based on the reduction of a phosphorus-molybdate complex with a reducing agent resulting in the subsequent formation of molybdenum blue. In 1972, Daly and Ertingshausen introduced a method for inorganic phosphorus which measured the unreduced phosphorus-molybdate complex. This procedure is a modification of the Daly and Ertingshausen technique.<sup>(2)</sup>

### TEST PRINCIPLE

Inorganic phosphorus reacts with ammonium molybdate in the presence of sulfuric acid to produce an unreduced phosphomolybdate complex. The increase in absorbance at 340 nm is directly proportional to the concentration of inorganic phosphorus in the sample.

### REAGENTS

Phosphorus Blank Reagent (R1): An acidic solution (pH < 1 at 25°C) containing a surfactant.

Phosphorus Molybdate Reagent (R2): An acidic solution (pH < 1 at 25°C) containing 2.8 mmol/L ammonium molybdate.

Phosphorus Calibrator: 1 x 15 mL of a solution containing 6.2 mg/dL (2 mmol/L) phosphorus. (For use with manual assays only.)

### WARNINGS & PRECAUTIONS FOR USE

#### Danger

Contains: Hydrochloric acid, Sulfuric acid

H314 – Causes severe skin burns and eye damage.

Prevention – P280 – Wear protective gloves/protective clothing/eye protection/face protection.

Do not breathe vapor/spray. Wash thoroughly after handling.

Response – P303 + P361 + P353 – IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 – Immediately call a POISON CENTRE or doctor/physician.

P301 + P330 + P331 – If swallowed: Rinse mouth. Do NOT induce vomiting.

P363 – Wash contaminated clothing before reuse.

P304 + P340 – If inhaled: Remove person to fresh air and keep comfortable for breathing.

Storage – P405 – Store locked up.

Disposal – P501 – Dispose of contents/container in accordance with local/regional/national/international regulations.

Avoid contact with skin and eyes.

See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are provided in a ready to use format.

Supplied reagent is stable at 18-26°C until expiry date stated on the labels. Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum. Serum should be separated promptly from the red blood cells because erythrocytes contain phosphate concentrations several times greater than those found in sera.<sup>(3)</sup>

### SAMPLE STORAGE

Specimens may be stored up to seven days at 2-8°C and up to six months when frozen at -20°C.

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(4)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, and hemolysis were evaluated for this phosphorus method on a Roche/Hitachi® 704 analyzer using a significance criterion of > 10% variance from control.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
mg/dL	mmol/L			
3.2	1.04	Hemoglobin	400 mg/dL	62.0 µmol/L
3.4	1.10	Bilirubin	40 mg/dL	684 µmol/L
3.2	1.03	Intralipid	600 mg/dL	1800 mg/dL (20 mmol/L) Simulated Triglycerides

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(5)</sup>

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Sekisui Diagnostics' Phosphorus-SL reagent.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material
3. Quality Control materials.

#### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:50 and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

## QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

## CALCULATIONS

The analyzer automatically calculates the inorganic phosphorus concentration of each sample.

## TEST LIMITATIONS

A sample with an inorganic phosphorus concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

## REFERENCE INTERVALS<sup>(1)</sup>

2.7-4.5 mg/dL (0.87-1.45 mmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 704 analyzer unless otherwise stated.

## RESULTS

Inorganic phosphorus concentration is reported as mg/dL (mmol/L).

## REPORTABLE RANGE (CLSI EP6)<sup>(4)</sup>

The linearity of the procedure described is 20 mg/dL (6.46 mmol/L). The lower limit of detection of the procedure described is 0.1 mg/dL (0.03 mmol/L). This data results in a reportable range of 0.1 to 20 mg/dL (0.03 to 6.46 mmol/L).

## PRECISION STUDIES (CLSI EP5)<sup>(4)</sup>

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV%	Concentration		Within Run SD		Within Run CV%
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	mg/dL	mmol/L	
3.3	1.07	0.06	0.02	1.7	3.3	1.07	0.05	0.02	1.5
7.5	2.42	0.10	0.03	1.4	7.7	2.49	0.04	0.01	0.6

Within run precision data was collected on two concentrations of control sera each run 20 times in a single assay.

## ACCURACY (CLSI EP9)<sup>(4)</sup>

The performance of this method (y) was compared with the performance of a similar phosphorus method (x) on a Roche/Hitachi® 704 analyzer. Forty-four patient serum samples ranging from 2.2-20.1 mg/dL (0.71-6.49 mmol/L) were tested and gave a correlation coefficient of 0.9998. Linear regression analysis gave the following equation:

$$\text{This method} = 0.9926 (\text{reference method}) + 0.00 \text{ mg/dL} (0.0 \text{ mmol/L}).$$

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**SEKISUI**  
DIAGNOSTICS

### The Americas

Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265

Fax: 902-628-6504

Email: questions@sekisuidiagnostics.com

peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

### International

Sekisui Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ, UK

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info@sekisuidiagnostics.com

**REFERENCES/ BIBLIOGRAFÍA**

1. Burtis, C.A., Ashwood, E.R., *Tietz Textbook of Clinical Chemistry, 2<sup>nd</sup> ed.* W.B. Saunders Co., Toronto, pp 1950-1952, p. 2202 (1994).
2. Daly, J.A. and Ertingshausen, G., *Direct Method for Determining Inorganic Phosphorus in Serum with the Centrifichem*, Clin. Chem. 18:263, (1972).
3. Kaplan, L., Pesce, A.J., *Clinical Chemistry, Theory, Analysis, Correlation, Third Edition*, Mosby, p. 552 (1996).
4. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Third Edition, Washington, 1990.

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KENT, ME16 0LQ, UK

United Kingdom/ Reino Unido

The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Fabricante

Consult instructions for use

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*In vitro* diagnostic medical device

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