UNSATURATED IRON BINDING CAPACITY (UIBC) ASSAY

CATALOGUE NUMBER: 153-01 **SIZE:** R1: 1 x 1000 mL
 153-10 R1: 1 x 100 mL + R2: 1 x 25 mL
 153-30 R1: 3 x 100 mL + R2: 1 x 75 mL
 153-50 R1: 1 x 500 mL + R2: 1 x 300 mL
 153-90 R1: 1 x 1000 mL + R2: 1 x 240 mL

INTENDED USE

For the IN VITRO quantitative measurement of unsaturated iron binding capacity (UIBC) in serum.

TEST SUMMARY

The measurement of unsaturated iron binding capacity (UIBC) in combination with serum iron is a useful diagnostic tool in the determination of various iron disorders. The combined value of UIBC and serum iron gives a value for the total iron binding capacity (TIBC). This represents the maximum concentration of iron that serum proteins can bind. Serum UIBC levels vary in disorders of iron metabolism, where iron binding capacities are often increased in iron deficiency and decreased in chronic inflammatory disorders or malignancies.⁽¹⁾

In 1970, Stookey⁽²⁾ reported the synthesis of 3-(2-pyridyl)-5,6-bis(4-phenylsulfonic acid)-1,2,4-triazine, monosodium salt (Ferrozine[®]) which complexed with ferrous iron to form a tris ferrozine/iron, Fe(Fz)₃ complex. The major advantages of ferrozine are the high molar absorptivity of the ferrous ferrozine complex (28,000), its water solubility, and stability over the pH range of 4-9. This assay uses a feroin type compound called 5,5'-(3-[2-pyridyl]-1,2,4-triazine-5,6 diyl)-bis-2-furansulfonic acid, disodium salt (Ferene[®]).^(3,4,5) This reagent is a superior iron chelating agent forming a Ferene[®] complex with ferrous iron with a maximum absorbance at 593 nm and a molar absorptivity of 35,500. The compound has a 27% higher molar absorption than ferrozine, absorbs at a longer wavelength, and has the other advantages of ferrozine, namely its solubility and stability.

TEST PRINCIPLE

Fe⁺⁺ (known) + Transferrin → Transferrin (Fe⁺⁺) + Fe⁺⁺ (excess)

Fe⁺⁺ (excess) + 3 Ferene → Ferrous Ferene (blue complex)

A known ferrous ion concentration incubated with serum binds specifically with transferrin at unsaturated iron binding sites.

Remaining unbound ferrous ions are measured with the ferene reaction. The difference between the amount of unbound iron and the total amount added to the serum is equivalent to the quantity bound to transferrin. This is the UIBC of the sample.

REAGENTS

Binding Buffer Reagent (R1): A buffer (pH 8.7 at 25°C) containing 12.8 µmol/L Ferrous Ammonium Sulfate, surfactants, stabilizers, and a preservative.

Color Reagent (R2): A solution containing 240 mmol/L Ascorbic Acid, 6.1 mmol/L Ferene, and a stabilizer.

WARNINGS & PRECAUTIONS FOR USE

European/International symbol required **Xn (Harmful)**.

S24/25: Avoid contact with skin and eyes.

R22: Harmful if swallowed.

R43: May cause sensitization by skin contact.

See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use.

Supplied reagent is stable until expiry date when stored at 2-8°C and protected from light. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum. The sample should be drawn in the morning following a 12 hour fast. Blood collecting glassware should be free of iron contamination.

SAMPLE STORAGE

If a sample is to be stored for longer than 8 hours, refrigeration at 2-8°C is recommended.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁶⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interference from icterus, lipemia, and hemoglobin was evaluated for this method on a Roche/Hitachi[®] 911 analyzer using a significance criterion of >10% variance from control.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
µg/dL	µmol/L			
323	57.6	Hemoglobin	400 mg/dL	62 µmol/L
236	42.2	Bilirubin	40 mg/dL	684 µmol/L
310	55.5	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

Copper is the only cation of the trace metals usually present in serum to form a colored complex with ferene. Copper interference with ferene is similar to that encountered with ferrozine and studied by Duffy and Gaudin.⁽⁷⁾ Ninety five percent of the copper interference is eliminated by chelation of free copper.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁸⁾

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' UIBC reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at the appropriate wavelength as per the instrument application.
2. Calibration material.
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:18.5, and wavelength readings of (primary/secondary) 600/700 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the UIBC concentration of each sample.

TEST LIMITATIONS

A sample with a UIBC concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS

UIBC reference intervals were determined from an apparently healthy population of 203 people following CLSI C28⁽⁶⁾ guidelines. Serum samples were analyzed on automated instrumentation at 37°C.

126-382 µg/dL (22.6-68.4 µmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 704 analyzer unless otherwise stated.

RESULTS

UIBC concentration is reported as µg/dL (µmol/L).

REPORTABLE RANGE (CLSI EP6)⁽⁶⁾

The linearity of the procedure described is 600 µg/dL (107.4 µmol/L). The lower limit of detection of the procedure described is 17 µg/dL (3.0 µmol/L). This data results in a reportable range of 17 to 600 µg/dL (3.0 to 107.4 µmol/L).

PRECISION STUDIES (CLSI EP5)⁽⁶⁾

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV (%)	Concentration		Within SD		Within Run CV (%)
µg/dL	µmol/L	µg/dL	µmol/L		µg/dL	µmol/L	µg/dL	µmol/L	
126	22.6	7.2	1.29	5.7	129	23.0	6.4	1.15	5.0
209	37.4	8.6	1.54	4.1	203	36.3	5.5	0.98	2.7

Within run precision data was collected on two concentrations of control sera each run 20 times in a single assay.

ACCURACY (CLSI EP9)⁽⁶⁾

The performance of this method (y) was compared with the performance of a similar UIBC method (x) on a Roche/Hitachi® 704 analyzer. Forty patient serum samples ranging from 94-437 µg/dL (16.8-78.1 µmol/L) were tested and gave a correlation coefficient of 0.9937. Linear regression analysis gave the following equation:

$$\text{This method} = 1.065 (\text{reference method}) + 2.0 \mu\text{g/dL} (0.4 \mu\text{mol/L}).$$

TRADEMARKS

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265

Fax: 902-628-6504

Email: questions@sekisuidiagnostics.com

peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

International
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK

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Fe⁺⁺ (exc. dent) + 3 Ferene → Ferene ferreux (complexe bleu)

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1. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., p. 1578 (1986).
2. Stookey, L.L., Ferrozine-A New Spectrophotometric Reagent for Iron, Anal. Chem. 42, 779 (1970).
3. Artiss, J.D., Vinogradov, S., Zak, B., Spectrophotometric Study of Several Sensitive Reagents for Serum Iron, Clin. Biochem. 14, 31-315 (1981).
4. Higgins, T., Novel Chromogen for Iron Determinations, Clin. Chem. 27, 1679 (1981).
5. Artiss, J.D., Strandbergh, D.R., Zak, B., Study of Continuous Flow Automation for Serum Iron on Comparing Several Sensitive Reagents, Microchemical Journal, 28, 275-284 (1983).
6. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
7. Duffy, J.R., Gaudin, J., Copper Interference in the Determination of Iron in Serum using Ferrozine, Clin. Biochem. 10, 122-123 (1977).
8. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Third Edition, Washington, 1990.

Authorized Representative/Représentant autorisé/Representante autorizado/
Rappresentante autorizzato/Autorisierte Vertreter/Geautoriseerde
Vertegenwoordiger/Yetkili Temsilci/Autoryzowany Przedstawiciel/
Εξουσιοδοτημένος αντιπρόσωπος

Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF
United Kingdom/Royaume-Uni/Reino Unido/Regno Unito/Vereinigtes
Königreich/Verenigd Koninkrijk/Ingiltere/Wielka Brytania/Ηνωμένο Βασίλειο
Tel (+44)(0)1732-220022
Fax (+44)(0)1732-220024

IN15310-16
June 7, 2011