MAGNESIUM REAGENT

CATALOGUE NUMBER: 175-12 **SIZE:** 4 x 125 mL

INTENDED USE

The Sekisui Magnesium Assay is intended for the quantitative measurement of magnesium in human serum and plasma (Lithium Heparin) on automated clinical chemistry analyzers. Magnesium measurements are used in the diagnosis and treatment of hypomagnesemia and hypermagnesemia. This device is intended for professional use and IN VITRO diagnostic use only.

TEST SUMMARY

Magnesium measurements are clinically important. Magnesium, along with potassium, is the most abundant intracellular cation. Magnesium is a required coenzyme in the metabolism of carbohydrates, lipids, and proteins.

Decreased magnesium concentrations, hypomagnesemia, may cause tremors, muscle irritability, increased blood pressure and heart rate.

Elevated magnesium concentrations, hypermagnesemia, may lead to impaired respiration, coma, and cardiac arrest if left uncorrected⁽¹⁾.

The method of choice for magnesium measurement is atomic absorption spectrophotometry⁽²⁾. This method is not practical for clinical laboratories so dye binding methods using methylthymol blue, titan yellow, calmagite, and xylidyl blue-1 are used instead. This method employs xylidyl blue-1 which allows the use of a single reagent which is rapid and precise.

TEST PRINCIPLE

Xylidyl blue-1 + Mg⁺⁺ → Mg-xylidyl blue complex (red)

Under the alkaline conditions of the Magnesium reagent, magnesium present in the sample, when mixed with the reagent, forms a red complex with the diazonium salt of xylidyl blue. The concentration of magnesium in the sample is directly proportional to the amount of xylidyl blue – magnesium complex formed and can be measured spectrophotometrically by the decrease in absorbance measured at 660 nm.

REAGENTS

Magnesium Reagent: A solution containing buffer (pH 11.2 at 25°C), 0.14 mmol/L xylidyl blue-1, 0.1 mmol/L EGTA, and a surfactant.

WARNINGS & PRECAUTIONS FOR USE

Contains: Ethanolamine

Danger

Hazard statements

Causes severe skin burns and eye damage. Causes serious eye damage. May cause respiratory irritation.

Precautionary statements

Prevention

Do not breathe mist or vapour. Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection.

Response

IF SWALLOWED: rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Remove immediately all contaminated clothing. Rinse skin with water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or physician.

See Material Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

Supplied reagent is stable at 2-8°C until the expiry date stated on the label. Stability claims are based on accelerated stability studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum and plasma. Lithium Heparin specimens were tested and found to be acceptable. Use only suitable tubes for collection and storage. Various manufacturers of sample collection systems may contain differing materials that may affect results. Follow all instructions provided by the collection system manufacturer.

SAMPLE STORAGE⁽³⁾

Magnesium is stable in serum for 7 days at 4-8°C and for 12 months at -20°C.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁴⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

EGTA prevents calcium interference by preferential chelation of the calcium.

Interferences from hemolysis, icterus, ascorbic acid and lipemia were evaluated for this magnesium method on a Roche/Hitachi® analyzer using a significance criterion of >10% variance from control. Interference data was collected in serum. Plasma is expected to be similar.

Concentration of Analyte		Substance Tested	Concentration of Interferent where Interference is Insignificant	
mg/dL	mmol/L			
2.1	0.86	Hemoglobin	1000 mg/dL	155 µmol/L
4.2	1.73		1000 mg/dL	155 µmol/L
7.9	3.25		1000 mg/dL	155 µmol/L
2.0	0.82	Conjugated Bilirubin	40 mg/dL	684 µmol/L
4.3	1.77		40 mg/dL	684 µmol/L
6.8	2.79	Unconjugated Bilirubin	40 mg/dL	684 µmol/L
2.0	0.82		40 mg/dL	684 µmol/L
4.2	1.73		40 mg/dL	684 µmol/L
6.8	2.79	Ascorbic Acid	3000 µg/dL	170 µmol/L
2.3	0.95		3000 µg/dL	170 µmol/L
4.4	1.81		3000 µg/dL	170 µmol/L
6.8	2.79	Intralipid	800 mg/dL	2400 mg/dL (27.1 mmol/L) Simulated Triglycerides
1.9	0.78		1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides
4.1	1.69		1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides
6.9	2.84		1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁵⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' Magnesium reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:117, and wavelength readings of 660/700 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the magnesium concentration of each sample.

TEST LIMITATIONS

A sample with a magnesium concentration exceeding the linearity limit should be diluted 1:1 with 0.9% saline and reassayed, incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽¹⁾

1.6-2.6 mg/dL (0.65-1.05 mmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 717 analyzer unless otherwise stated.

RESULTS

Magnesium concentration is reported as mg/dL (mmol/L).

To convert magnesium results to mEq/L, use the following conversion factor:⁽⁶⁾
 $\text{mmol/L} \times 2 = \text{mEq/L}$ ($\text{mEq/L} \times 0.5 = \text{mmol/L}$)

REPORTABLE RANGE (CLSI EP6)⁽⁴⁾

The linearity of the procedure described is 8.0 mg/dL (3.29 mmol/L). The limit of blank and limit of detection of the procedure described are 0.2 mg/dL (0.08 mmol/L) and 0.3 mg/dL (0.12 mmol/L), respectively. The limit of quantitation of the procedure described is 0.3 mg/dL (0.12 mmol/L). This data results in a reportable range of 0.3 mg/dL (0.12 mmol/L) to 8.0 mg/dL (3.29 mmol/L). The reportable range will depend on the sample to reagent ratio used.

PRECISION STUDIES (CLSI EP5)⁽⁴⁾

Total precision data was collected on three concentrations of serum in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Concentration		Within Run SD		With-in Run CV %
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	mg/dL	mmol/L	
0.8	0.33	0.05	0.021	6.1	0.9	0.37	0.06	0.024	7.0
2.0	0.82	0.07	0.029	3.7	2.0	0.82	0.04	0.016	2.2
4.7	1.93	0.15	0.062	3.2	4.9	2.01	0.08	0.033	1.6

Within run precision data was collected on three concentrations of serum each run 20 times in a single assay.

ACCURACY (CLSI EP9)⁽⁴⁾

The performance of this method (y) was compared with the performance of a similar xylidyl blue method (x) on a Roche/Hitachi® 717 analyzer. One hundred patient serum samples ranging from 0.5 mg/dL to 7.3 mg/dL (0.21-2.92 mmol/L) gave a correlation coefficient of 0.985. Linear regression analysis gave the following equation:
This method = 0.981 (reference method) - 0.0 mg/dL (0.02 mmol/L).

The performance of this method (plasma) was compared with the performance of this method (serum) on a Roche/Hitachi® 717 analyzer. Fifty-four (54) patient serum and plasma samples ranging from 0.5 mg/dL to 5.4 mg/dL (0.21-2.21 mmol/L) gave a correlation coefficient of 0.9903. Linear regression analysis gave the following equation:
This method (plasma) = 1.009 (serum) - 0.1 mg/dL (0.05 mmol/L).

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

SEKISUI
DIAGNOSTICS

The Americas
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
peidiagnosticttechnical@sekisuidiagnostics.com

International
Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, UK
Email: info@sekisuidiagnostics.com

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peidiagnosticttechnical@sekisuidiagnostics.com

Internacional

Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT, ME19 4AF, RU

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info@sekisuidiagnostics.com

www.sekisuidiagnostics.com

**Definitions for Symbols/ Définitions des Symboles/
Definición de los Símbolos**

This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Manufacturer
Fabricant
Fabricante

Consult instructions for use
Consulter les directives d'utilisation
Consulte las instrucciones de uso

In vitro diagnostic medical device
Appareil médical de diagnostic *in vitro*.
Dispositivo médico para el diagnóstico *in vitro*

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Temperature limitation
Limite de température
Límites de temperatura

REFERENCES/RÉFÉRENCES/ REFERENCIAS

1. Tietz, N.W., (Editor) Clinical Guide to Laboratory Tests, W.B. Saunders Company, Philadelphia (1983) p.338.
2. Ratge, D., Kohse, K.P., and Wisser, H.; Measurement of Magnesium in Serum and Urine with a Random Access Analyzer by use of a Modified Xylidyl Blue-1 Procedure, Clinical Chimica Acta, 159 (1986) 197-203.
3. World Health Organization; Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/Lab/99.1 Rev. 2, 2002.
4. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, Third Edition, 1990.
6. Tietz, N.W. (Editor) Clinical Guide to Laboratory Tests, W.B. Saunders Company, Philadelphia (1995) p. 418.

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Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF
United Kingdom/ Reino Unido
Tel (+44)(0)1732-220022
Fax (+44)(0)1732-220024

IN17512-0
January 2, 2012