EN

CARBON DIOXIDE L3K[®] ASSAY


CATALOGUE NUMBER: 299-30 **SIZE:** 6 x 30 mL
 299-17 5 x 100 mL
 299-50 4 x 500 mL
 299-55 1 x 500 mL

INTENDED USE

For the IN VITRO quantitative measurement of carbon dioxide in serum and plasma.

TEST SUMMARY


Diagnostic determination of serum or plasma carbon dioxide (CO₂) is used in conjunction with other clinical and laboratory testing for the evaluation of acid-base balance. Blood pH is buffered by the presence of the weak carbonic acid (H₂CO₃) and conjugate base bicarbonate (HCO₃⁻). As CO₂ concentration in blood increases, so does hydrogen ion (H⁺) concentration. Respiration rate, which is controlled by CO₂ concentration in the blood, compensates by increasing or decreasing in an attempt to maintain equilibrium by promoting CO₂ elimination or retention respectively.

High CO₂ content is often indicative of respiratory acidosis or metabolic alkalosis. Low CO₂ content is often indicative of respiratory alkalosis or metabolic acidosis. Further laboratory testing is necessary to differentiate between metabolic and respiratory conditions.⁽¹⁾


Classic techniques for the measurement of carbon dioxide (CO₂) involve the addition of acid to liberate the carbon dioxide and the measurement of carbon dioxide thus released by manometric, volumetric, or titrimetric techniques. These procedures are both time consuming and cumbersome. Sekisui Diagnostics' Carbon Dioxide L3K[®] assay is an enzymatic procedure, employing phosphoenolpyruvate carboxylase (PEPC)⁽²⁾ and a stabilized NADH analog,⁽³⁾ which is easy to use and applicable to routine laboratory instrumentation.

TEST PRINCIPLE

PEPC catalyses the first reaction which produces oxaloacetate. In the presence of MDH, the reduced cofactor is oxidized by oxaloacetate. The decrease in concentration of the reduced cofactor is monitored between 405 and 415 nm and is proportional to the total carbon dioxide concentration in the sample.

PEPC is specific for the bicarbonate ion (HCO₃⁻) and its action disturbs the following equilibrium which results in conversion of the CO₂ to HCO₃⁻.

REAGENTS

Carbon Dioxide Reagent: A solution containing buffer (pH 7.5 at 25°C), 12.5 mmol/L PEP, >400 U/L PEPC (microbial), >4100 U/L malate dehydrogenase (mammalian), 0.6 mmol/L NADH analog, activators, stabilizers, a surfactant, and a preservative.

WARNINGS & PRECAUTIONS FOR USE

IVD

For *In Vitro* Diagnostic Use.

Rx ONLY

When using, do not eat or drink.
 Avoid contact with skin and eyes.

Warning: Suspected of damaging fertility or the unborn child. May cause damage to organs (central nervous system, kidney) through prolonged or repeated exposure. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Do not breathe mist/vapors/spray. Wear protective gloves/protective clothing/eye protection/face protection. If exposed or concerned: Get medical advice/attention.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use.
 Supplied reagent stable at 2-8°C until expiry date.
 Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum or lithium heparinized plasma. The specimen should be promptly separated from the clot or cells and stored tightly sealed to prevent the loss of carbon dioxide.⁽⁶⁾

SAMPLE STORAGE

Samples should be stored at 2-8°C.⁽⁷⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁴⁾

Cross Contamination studies have not been performed on automated instruments. Certain reagent / instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia and hemolysis were evaluated for this carbon dioxide method on a Roche/Hitachi[®] 717 analyzer using a significance criterion of >10% variance from control. Data for plasma is similar to that presented for serum.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
26.3 mmol/L	26.3 mEq/L	Hemoglobin	600 mg/dL	93 µmol/L
27.1 mmol/L	27.1 mEq/L	Bilirubin	40 mg/dL	684 µmol/L
24.6 mmol/L	24.6 mEq/L	Intralipid	1000 mg/dL	3000 mg/dL (34 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostic's studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁵⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' Carbon Dioxide L3K[®] Reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a fixed rate test mode, with a sample to reagent ratio of 1:100, and a wavelength reading of 405 nm or 415 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the carbon dioxide concentration of each sample.

TEST LIMITATIONS

A sample with a carbon dioxide concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS

22-29 mmol/L (22-29 mEq/L)⁽⁶⁾

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on Roche/Hitachi® 717 analyzer unless otherwise stated.

RESULTS

Carbon dioxide concentration is reported as mmol/L (mEq/L).

REPORTABLE RANGE (CLSI EP6)⁽⁴⁾

The linearity of the procedure described is 50.0 mmol/L (mEq/L). The lower limit of detection of the procedure described is 2.9 mmol/L (mEq/L). This data results in a reportable range of 2.9 to 50.0 mmol/L (mEq/L).

PRECISION STUDIES (CLSI EP5)⁽⁴⁾

Total precision data was collected on two concentrations of control sera in triplicate over five runs, for a total of 15 sample values.

Concentration mmol/L (mEq/L)	Total SD mmol/L (mEq/L)	Total CV%	Concentration mmol/L (mEq/L)	Within Run SD mmol/L (mEq/L)	Within Run CV%
14	0.24	1.7	14	0.18	1.3
25	0.27	1.1	24	0.40	1.7

Within run precision data was collected on two concentrations of control sera, each run 20 times in a single assay.

ACCURACY (CLSI EP9)⁽⁴⁾

The performance of this method (y) was compared with the performance of a similar carbon dioxide method (x) on a Cobas Mira® analyzer. Forty-five serum samples ranging from 15.4 to 43.7 mmol/L (mEq/L) were tested and gave a correlation coefficient of 0.9861. Linear regression analysis gave the following equation:

$$\text{This method} = 1.00(\text{reference method}) - 0.42 \text{ mmol/L (mEq/L)}$$

The performance of this method (y) with plasma samples (y) was compared to this method with serum samples (x) on an Advia® 1650 analyzer.

Twenty-five patient serum and plasma samples ranging from 17.1 to 25.3 mmol/L (mEq/L) gave correlation coefficient of 0.9636. Linear regression analysis gave the following equation:

$$\text{This method (plasma)} = 0.950[\text{this method(serum)}] + 1.35 \text{ mmol/L (mEq/L)}$$

TRADEMARKS

L3K is a registered trademark of Sekisui. All other trademarks, brands, product names are the property of their respective companies.

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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1. Tietz, N.W. (Ed.), *Fundamentals of Clinical Chemistry*, W.B. Saunders Co., Toronto, 625-636, 937 (1970).
2. Norris, K.A., Atkinson, A.R., Smith, W.G., Clin. Chem. 21 (1975).
3. US Patent #5,801,006
4. *CLSI Guidelines and Standards*, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Third Edition, Washington (1990).
6. Tietz, N.W. (Editor), *Textbook of Clinical Chemistry*, W.B. Saunders Co., Philadelphia (1986).
7. World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations. Geneva: World Health Organization; 2002:39.

PATENT NUMBER

US Patent No. 5,801,006

EC REP

MDSS GmbH
Schiffgraben 41,
30175 Hannover, Germany
Phone: (+49)-511-6262 8630
Fax: (+49)-511-6262 8633

The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.

IN29930-18
September 25, 2019