EN

## MAGNESIUM REAGENT

**CATALOGUE NUMBER:** 175-12      **SIZE:** 4 x 125 mL

### INTENDED USE

The Sekisui Magnesium Assay is intended for the quantitative measurement of magnesium in human serum and plasma (Lithium Heparin) on automated clinical chemistry analyzers. Magnesium measurements are used in the diagnosis and treatment of hypomagnesemia and hypermagnesemia. This device is intended for professional use and IN VITRO diagnostic use only.

### TEST SUMMARY

Magnesium measurements are clinically important. Magnesium, along with potassium, is the most abundant intracellular cation. Magnesium is a required coenzyme in the metabolism of carbohydrates, lipids, and proteins.

Decreased magnesium concentrations, hypomagnesemia, may cause tremors, muscle irritability, increased blood pressure and heart rate.

Elevated magnesium concentrations, hypermagnesemia, may lead to impaired respiration, coma, and cardiac arrest if left uncorrected<sup>(1)</sup>.

The method of choice for magnesium measurement is atomic absorption spectrophotometry<sup>(2)</sup>. This method is not practical for clinical laboratories so dye binding methods using methylthymol blue, titan yellow, calmagite, and xylydyl blue-1 are used instead. This method employs xylydyl blue-1 which allows the use of a single reagent which is rapid and precise.

### TEST PRINCIPLE

$\text{Xylydyl blue-1} + \text{Mg}^{++} \rightarrow \text{Mg-xylydyl blue complex (red)}$

Under the alkaline conditions of the Magnesium reagent, magnesium present in the sample, when mixed with the reagent, forms a red complex with the diazonium salt of xylydyl blue. The concentration of magnesium in the sample is directly proportional to the amount of xylydyl blue – magnesium complex formed and can be measured spectrophotometrically by the decrease in absorbance measured at 660 nm.

### REAGENTS

Magnesium Reagent: A solution containing buffer (pH 11.35 at 25°C), ~0.14 mmol/L xylydyl blue-1, 0.1 mmol/L EGTA, and a surfactant.

### WARNINGS & PRECAUTIONS FOR USE

**IVD**

For *In Vitro* Diagnostic Use.

**Rx ONLY**

#### Danger

Contains: Ethanolamine  
H314 - Causes severe skin burns and eye damage.  
H335 – May cause respiratory irritation.  
Precautionary statements  
Prevention P260 – Do not breathe the mist or vapour.  
P264 – Wash thoroughly after handling.  
P280 – Wear protective gloves/protective clothing/eye protection/face protection.  
Response – P301 + P330 + P331 – IF SWALLOWED: rinse mouth. Do NOT induce vomiting.  
P303 + P361 + P353 – IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.  
P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P304 + P340 – IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing.  
P310 – Immediately call a POISON CENTRE or doctor/physician.  
Storage – P405 – Store locked up.  
Disposal – P501 – Dispose of contents/container in accordance with local/regional/national/international regulations.

See Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

Supplied reagent is stable at 2-8°C until the expiry date stated on the label. Stability claims are based on accelerated stability studies.

### REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum and plasma. Lithium Heparin specimens were tested and found to be acceptable. Use only suitable tubes for collection and storage. Various manufacturers of sample collection systems may contain differing materials that may affect results. Follow all instructions provided by the collection system manufacturer.

### SAMPLE STORAGE<sup>(3)</sup>

Magnesium is stable in serum for 7 days at 4-8°C and for 12 months at -20°C.

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(4)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

EGTA prevents calcium interference by preferential chelation of the calcium.

Interferences from hemolysis, icterus, ascorbic acid and lipemia were evaluated for this magnesium method on a Roche/Hitachi® analyzer using a significance criterion of >10% variance from control. Interference data was collected in serum. Plasma is expected to be similar.

Concentration of Analyte		Substance Tested	Concentration of Interferent where Interference is Insignificant	
mg/dL	mmol/L			
2.1	0.86	Hemoglobin	1000 mg/dL	155 µmol/L
4.2	1.73		1000 mg/dL	155 µmol/L
7.9	3.25		1000 mg/dL	155 µmol/L
2.0	0.82	Conjugated Bilirubin	40 mg/dL	474 µmol/L
4.3	1.77		40 mg/dL	474 µmol/L
6.8	2.79		40 mg/dL	474 µmol/L
2.0	0.82	Unconjugated Bilirubin	40 mg/dL	684 µmol/L
4.2	1.73		40 mg/dL	684 µmol/L
6.8	2.79		40 mg/dL	684 µmol/L
2.3	0.95	Ascorbic Acid	3000 µg/dL	170 µmol/L
4.4	1.81		3000 µg/dL	170 µmol/L
6.8	2.79		3000 µg/dL	170 µmol/L
1.9	0.78	Intralipid	800 mg/dL	2400 mg/dL (27.1 mmol/L) Simulated Triglycerides
4.1	1.69		1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides
6.9	2.84		1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(5)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Sekisui Diagnostics' Magnesium reagent.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

#### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:117, and wavelength readings of 660/700 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

## CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

## QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

## CALCULATIONS

The analyzer automatically calculates the magnesium concentration of each sample.

## TEST LIMITATIONS

A sample with a magnesium concentration exceeding the linearity limit should be diluted 1:1 with 0.9% saline and reanalyzed, incorporating the dilution factor in the calculation of the value.

## REFERENCE INTERVALS<sup>(1)</sup>

1.6-2.6 mg/dL (0.66-1.07 mmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 717 analyzer unless otherwise stated.

## RESULTS

Magnesium concentration is reported as mg/dL (mmol/L).

To convert magnesium results to mEq/L, use the following conversion factor:<sup>(6)</sup>  
 $\text{mmol/L} \times 2 = \text{mEq/L}$  ( $\text{mEq/L} \times 0.5 = \text{mmol/L}$ )

## REPORTABLE RANGE (CLSI EP6)<sup>(4)</sup>

The linearity of the procedure described is 8.0 mg/dL (3.29 mmol/L). The limit of blank and limit of detection of the procedure described are 0.2 mg/dL (0.08 mmol/L) and 0.3 mg/dL (0.12 mmol/L), respectively. The limit of quantitation of the procedure described is 0.3 mg/dL (0.12 mmol/L). This data results in a reportable range of 0.3 mg/dL (0.12 mmol/L) to 8.0 mg/dL (3.29 mmol/L). The reportable range will depend on the sample to reagent ratio used.

## PRECISION STUDIES (CLSI EP5)<sup>(4)</sup>

Total precision data was collected on three concentrations of serum in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Concentration		Within Run SD		With-in Run CV %
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	mg/dL	mmol/L	
0.8	0.33	0.05	0.021	6.1	0.9	0.37	0.06	0.024	7.0
2.0	0.82	0.07	0.029	3.7	2.0	0.82	0.04	0.016	2.2
4.7	1.93	0.15	0.062	3.2	4.9	2.01	0.08	0.033	1.6

Within run precision data was collected on three concentrations of serum each run 20 times in a single assay.

## ACCURACY (CLSI EP9)<sup>(4)</sup>

The performance of this method (y) was compared with the performance of a similar xylidyl blue method (x) on a Roche/Hitachi® 717 analyzer. One hundred patient serum samples ranging from 0.5 mg/dL to 7.3 mg/dL (0.21-3.00 mmol/L) gave a correlation coefficient of 0.985. Linear regression analysis gave the following equation:

$$\text{This method} = 0.981 (\text{reference method}) - 0.0 \text{ mg/dL} (0.02 \text{ mmol/L}).$$

The performance of this method (plasma) was compared with the performance of this method (serum) on a Roche/Hitachi® 717 analyzer. Fifty-four (54) patient serum and plasma samples ranging from 0.5 mg/dL to 5.4 mg/dL (0.21-2.22 mmol/L) gave a correlation coefficient of 0.9903. Linear regression analysis gave the following equation:

$$\text{This method (plasma)} = 1.009 (\text{serum}) - 0.1 \text{ mg/dL} (0.04 \text{ mmol/L}).$$

The information presented above is based on results from Sekisui Diagnostics studies and is current at the date of publication.

All trademarks, brands, product names and trade names are the property of their respective companies.

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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Manufacturer  
Fabricant  
Fabricante

Consult instructions for use  
Consulter les directives d'utilisation  
Consulte las instrucciones de uso

*In vitro* diagnostic medical device  
Appareil médical de diagnostic *in vitro*.  
Dispositivo médico para el diagnóstico *in vitro*

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Temperature limitation  
Limite de température  
Límites de temperatura

For use by or on the order of a physician only (applicable to USA classification only)

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**REFERENCES/RÉFÉRENCES/ REFERENCIAS**

1. Tietz, N.W., (Editor) Clinical Guide to Laboratory Tests, W.B. Saunders Company, Philadelphia (1983) p.338, 339.
2. Ratge, D., Kohse, K.P., and Wisser, H.; Measurement of Magnesium in Serum and Urine with a Random Access Analyzer by use of a Modified Xylidyl Blue-1 Procedure, Clinical Chimica Acta, 159 (1986) 197-203.
3. World Health Organization; Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/Lab/99.1 Rev. 2, 2002.
4. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, Third Edition, 1990.
6. Tietz, N.W. (Editor) Clinical Guide to Laboratory Tests, W.B. Saunders Company, Philadelphia (1995) p. 418.

MDSS GmbH  
Schiffgraben 41,  
30175 Hannover, Germany  
**Phone:** (+49)-511-6262 8630  
**Fax:** (+49)-511-6262 8633

The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.

IN17512-5  
September 25, 2019